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- WHO Influenza Reagent Kit for the Identification of Influenza A/H5 Virus (1997–98, 2003 or 2004 version). The reagents in this kit for the immunofluorescence assay include:
  - influenza type A/H5-specific monoclonal antibody pool
  - influenza A type-specific and influenza B type-specific monoclonal antibody pools
  - influenza A/H1 and an A/H3 subtype specific monoclonal antibodies
- Anti-mouse IgG FITC conjugate
- Microscope slides
- Cover slips, 24 x 60 mm
- Mountant
- Acetone
- Immunofluorescence microscope.

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- Madin-Darby Canine Kidney cells (MDCK). ATCC CCL34.
- WHO Influenza Reagent Kit for the Identification of Influenza A/H5 Virus.  
Reagents for identification of A/H5 in cell culture includes:
  - influenza A/H5 control antigen: inactivated virus
  - goat serum to A/Tern/South Africa/61/ H5
  - chicken pooled serum to A/Goose/Hong Kong/437-4/99
- WHO Influenza Reagent Kit (Annually distributed)  
A (H1N1) and A (H3N2) reference antigens and reference antisera.
- Receptor-destroying enzyme (RDE).
- Red blood cells ( chicken, turkey, human type O, or guinea-pig red blood cells) in Alsever's solution.

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「WHO laboratory biosafety guidelines for handling specimens suspected of  
containing highly pathogenic avian influenza A/H5 virus」

原文：[http://www.who.int/csr/disease/avian\\_influenza/guidelines/handlingspecimens/en/](http://www.who.int/csr/disease/avian_influenza/guidelines/handlingspecimens/en/)

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HA gene primers use by the Hong Kong Government Virus Unit modified from  
Yuen et al., 1998:

H5-1: GCC ATT CCA CAA CAT ACA CCC

H5-2: TAA ATT CTC TAT CCT CCT TTC CAA

Expected product size: 358bp

NA gene primers used by the Hong Kong Government Virus Unit modified from Wright et al 1995:

N1-1: TTG CTT GGT CGG CAA GTG C

N1-2: CCA GTC CAC CCA TTT GGA TCC

Expected product size: 615bp

- Positive control (Obtained upon request from a WHO Collaborating Centre for Influenza)
- 10X PCR buffer
- Adjustable pipettes, 10, 20 and 100  $\mu$ l
- Disposable filter tips
- Microcentrifuge, adjustable to 13 000 rpm
- Vortex mixer
- Thermocycler
- Agarose gel casting tray , electrophoresis chamber and power supply
- UV-light box or hand-held UV light (302 nm)

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a. Aliquot 45  $\mu$ l of master mix into each 0.2-ml PCR tube.

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Fouchier RA et al. (2000). Detection of influenza A viruses from different species by PCR amplification of conserved sequences in the matrix gene. *Journal of Clinical Microbiology*, 38:4096–4101.

Lennette EH, Schmidt NJ, eds (1979). *Diagnostic procedures for viral, rickettsial and chlamydial infections*, 5th ed. Washington, DC, American Public Health Association.

Nicholson KG, Wood JM, Zambon M (2003). Influenza. *Lancet*, 362:1733–1745.

WHO (2002). *WHO manual on animal influenza diagnosis and surveillance*. Geneva, World Health Organization (document WHO/CDS/CSR/NCS/2002.5, available at: <http://www.who.int/csr/resources/publications/influenza/en/whocdscsrncs2002.5.pdf>).

Wright KE et al. (1995). Typing and subtyping of influenza viruses in clinical samples by PCR. *Journal of Clinical Microbiology*, 33:1180–1184.

Yuen KY et al. (1998). Clinical features and rapid viral diagnosis of human disease associated with avian influenza A H5N1 virus. *Lancet*, 351:467–471.